

**YADU**  
M E D I C A L

# Surgical Drapes

*Establish aseptic protection barrier*

YADU surgical drapes series



HeNan YADU Industrial Co., Ltd.





## HeNan YADU Industrial Co., Ltd.

### Company Profile

Established in 1992, YADU is a dedicate company of developing, producing and suppling consumable medical devices, self-owned and operate production lines of nonwoven converting, gauze spinning, waving & sawing, injection moulding, extrusion and many more. We supply more than 100 categories, 10,000 SKU' s medical & healthcare products to over 8,000 OEM clients, distributors and hospitals.

**5**

Facilities in China & Bangladesh

**5,000**  
Employees

**700,000m<sup>2</sup>**  
Total Premises



## Surgical Drapes Workshop





## Product Superiority

- Boost efficiency with versatile options for every procedure.
- Promote safety with high-tech materials that help protect against infection.
- Convenient to locate the incision by designing scientific and reasonable fenestration for different operations
- Hand-applied premium tape and incise film offer precise coverage to help prevent pooling.
- Wide anesthesia screens accommodate larger patients comfortably.
- Drapes are not made with natural rubber latex unless indicated in the description.

## Surgical Drapes

Basic and universal procedures	Extremity procedures
Laparotomy procedures	Arthroscopy procedures
Cardiovascular procedures	Hip procedures
Cystoscopy procedures	General orthopedic procedures
Obstetrics procedures	Shoulder procedures
EENT/opthalmology procedures	Lithotomy/laparoscopy procedures

We produce surgical drapes for above surgical procedures. Surgical drapes shown on this brochure are only part of our products. For drapes not displayed or customized drapes, please discuss with our sales representative.



# Regulatory

Both design & construction of YADU plants and recruitment of our team are on account of international standards and regulatory requirements. Our operation are strictly based on QMS and GMP practices. We have been granted with certificates of CE, FDA 510k, BSCI, ISO 13485, ISO 9001 and many more.



Email: [info@yadu-med.com](mailto:info@yadu-med.com)


Web/English: <http://www.yadu-med.com>

Address: 66 Yangze Rd, Changyuan, Xinxiang, Henan, China




## Cesarean Drape

- Fenestration, full incise, 23 x 35.5cm
- Reinforcement below pouch, 43 x 56cm
- Armboard covers
- 360-degree fluid control pouch with contoured foam ring, Cord-Holding tab and ports
- Cord-Holding Tabs




## Arthroscopy Drape

- Heavyweight, premium fabric
- Absorbent reinforcement
- Dual fenestration, circular elastic, 6cm
- Anesthesia screen
- Armboard covers
- Cord-holding tabs
- Fluid control pouch with malleable wire drawstring pouch and port




## Angiography Drape

- Absorbent reinforcement, 73 x 112cm
- Fenestration, dual circular, 15cm with circular incise aperture, 9cm
- Single clear side panel on patient's right side






## Laparotomy Drape




- Absorbent reinforcement, 77.5 x 77.5 cm
- Fenestration, 10 x 30.5cm
- Instrument pad
- Armboard covers
- Cord-Holding Tabs

## TUR Drape



- Absorbent reinforcement, 51 x 48cm
- Fenestration with adhesive, 6 x 16.5cm
- Anesthesia screen
- Attached leggings
- Fluid control pouch with filtering screen, suction port, drawstring feature and side trough
- Cord-holding tabs

## Hip drape



- Absorbent reinforcement, 54 x 145cm
- Fenestration, pear-shaped elastic, 17 x 21cm
- Anesthesia screen
- Armboard covers
- Cord-Holding Tabs
- Clear leg positioning pouches, 138 x 76cm




## Under-buttocks drape

- with fluid control
- pouch, filtering screen, drainage port and
- reinforcement




## Split Drape

- Absorbent reinforcement
- Split with adhesive, 13 x 203cm
- Armboard covers



## Ophthalmic Drape

- Dual saddlebag fluid control pouches with malleable strips
- Malleable breather-bar
- Cord-holding tab
- Full incise, 10 x 13cm






## Side Drape

- Adhesive along reinforced edge



## Mayo Cover

- Helps ease the application of O.R. accessories
- High-quality lamination material provides exceptional barrier protection



## Leggings

- Our standard drape collection features a top-of-the-line SMS fabric for exceptionally high barrier strength
- Made of 100% polypropylene, SMS drapes are soft and flexible



## Stockinettes

- Impervious outer layer for protection; fibrous inner layer for comfort
- Strong barrier helps prevent fluid strikethrough
- Used for handling limbs in the OR, primarily during orthopedic procedures



## Back table cover

- Top-of-the-line SMS fabric is composed of 5 strong fiber layers for exceptionally high barrier strength
- Smaller fibers and a more uniform fiber layer enhance softness and drape

