

YADU
M E D I C A L

Email: info@yadu-med.com

Web/English: <http://www.yadu-med.com>

Address: 66 Yangze Rd, Changyuan, Xinxiang, Henan, China



YADU
M E D I C A L

Since 1992



YADU Industrial Co., Ltd.



建康中國 壯舉同心 守正出奇 行穩致遠

YADU 亚都集团

亚都集团 YADU GROUP



About us

YADU Medical is one of the largest and leading manufacturers of medical consumables. Established in 1992, YADU dedicate to developing and producing medical consumables for over 30 years. Our 5,000 employees runs laboratories, cleanrooms, environment controlled workshops, sterilization units and warehouses in five production facilities located in both China and Bangladesh, total land size is more than 70,000 square meters. Our products are supplied to more than 20 countries and over 8,000 local hospitals.

5

Factories

72,000
m²

Land size


5,000

Employees

Contents

01	Surgical Gown	07	Universal Pack
02	Isolation Gown		Head and Neck Pack
	Shoe cover	08	Laparotomy pack
	Cap		Laparoscopy pack
03	Top drape with armboard covers	09	Cesarean Pack
	Side drape		Delivery pack
	Bottom drape	10	Extremity Pack
04	Back Table Cover		Joint pack
	Fenestrated drape	11	Orthopedic pack
	Laparotomy T-drape with incise and instrument pad		Arthroscopy Pack
05	CSR Wrap	12	Angiography pack
	Underpad		TUR/Urology Pack
06	Adult Diaper	13	Dialysis Kit
	Wipes		Catheter Kit
		14	Debridement Kit
			Suture Kit
		15	Dressing Kit
			I.V. Kit
		16	First aid Kit
			First aid Kit
		17	Foley Catheter
		18	Connecting Catheter
		19	Oxygen Tube
		20	Breathing Circuit
		21	Vagunal Speculum
			Surgical Gloves
		22	Gauze swab
			Medical Face Mask
		23	Medical Particulate Respirator
			Coverall







Surgical Gowns

EN13795 Standard Performance & High Performance
 Individually Packed
 AAMI Level 2, Level 3 & Level 4

SPP+PE Reinforcement
 Blue, Green

Double Wrapping with SMS Sterile Field
 SMS, BVB or Spunlace
 Spunlace Towel, Paper Towel or Hulk Towel

S M L XL XXL



Isolation Gowns

Level 1, Level 2, Level 3
 SMS, SPP+PE, SPP
 Blue, Yellow

S M L XL XXL

Caps

Blue, White
 SMS, SPP+PE, SPP
 21", 24"




Shoe Covers

Blue, White, Green
 SMS, SPP+PE, SPP
 15.5", 16.5", 18"




Surgical Drapes




Top Drape with Armboard Covers

Absorbent reinforcement, 24"x16"(61x41cm)
Split, 24"x19"(61x48cm); with 2"(5cm)wide adhesive
Hook-and-loop line holders




Back Table Covers

SMMS, 2-layer Lamination
150x200cm or Customized




Side Drapes

Absorbent reinforcement, 8"x40"(20x102cm)
Adhesive along reinforced edge
Fabric trough, 4"x40"(10x102cm)
Clipping tabs




Fenestrated Drapes

Absorbent reinforcement, 12"x13"(30x33cm)
Fenestration, 4"x6"(10x15cm)



Bottom Drapes

Absorbent reinforcement, 29"x16"(74x41cm)
Adhesive along reinforced edge
Hook-and-loop line holders



Laparotomy T-drape with Incise and Instrument Pads

Absorbent reinforcement, 28"x23"(71x58cm)
Fenestration with incise, 9"x14"(23x36cm)
Instrument pad
Contoured armboard covers
Hook-and-loop line holders




CSR Wraps

- Strength and resistance to punctures, tears and abrasions
- Anti-static coating reduces static “shocks” to clinicians
- Low linting
- Soft and easy to fold
- Excellent drapability and low memory
- Excellent microbial barrier
- Effective sterilant penetration



Underpads

- Designed to lock wetness with light absorbency
- Fluff / polymer core with superior mat configuration.
- Heat-sealed, non-woven edges are gentle on skin.
- Great protection for general procedures and wound care.
- Soft, non-woven bonded topsheet promotes rapid fluid acquisition.



Adult Diapers

- The core with odor protection wicks away fluid to help promote dryness and maintain skin integrity
- Skin safe closures provide secure, repeated refastening
- Ultra soft cloth-like back sheet provides a discreet, garment-like feel and breathable side panels allow airflow for comfort
- Wetness indicator changes color when the garment has been soiled
- Moderate absorbency for total incontinence protection

- White or customized colour and/or printing
- S, L, 3XL, 5XL or Customized



Wipes

- Full sized, high weight, durable wet wipes with a touch of soothing aloe are ideal for everyday clean-ups and incontinence care
- The pH-balanced, hypoallergenic and alcohol-free wipes provide exceptional skin care
- Single patient use wipes help to reduce the risk of cross contamination
- Produce and supply both Wet Wipe and Dry Wipe

- Spunlace Nonwoven, EDI Pure Water, Chamomile Extract, Vitamin E 14x18cm, 14x20cm, 15x20cm, or Customized



Universal Packs

NO.	BOM	Qty	Unit
1	Mayo Stand Cover 80x145cm	1	EA
2	Op-tapes 10x50cm	2	EA
3	Hand Towels 40x40cm	2	EA
4	Side Drapes 75x90cm	2	EA
5	Top Drape 150x240cm	1	EA
6	Bottom Drape 150x180cm	1	EA
7	Wrapper 100x100cm	1	EA
8	Reinforced Surgical Gowns L	2	EA
9	Instrument Table Cover 150x200cm	1	EA



Laparotomy Packs

NO.	BOM	Qty	Unit
1	Instrument Table Cover 150x220cm	1	EA
2	Abdominal Drape 230x270cm	1	EA
3	Kidney Dish 700ml	1	EA
4	Mayo Stand Cover 80x145cm	1	EA
5	Side Drape 75x90cm	1	EA
6	Needle Counter 25 counts	1	EA
7	Hand Towels 30x40cm	4	EA
8	Standard Surgical Gowns L	2	EA



Head and Neck Packs

NO.	BOM	Qty	Unit
1	Hand Towels 30x34cm	4	EA
2	Standard Surgical Gown 125x150cm	1	EA
3	Reinforced Surgical Gown 125x150cm	1	EA
4	Split Drape 180x290cm	1	EA
5	Head Drape 100x100cm	1	EA
6	Side Drape 130x115cm	1	EA
7	Suture Bag 16x29cm	1	EA
8	Instrument Table Cover 150x190cm	1	EA
9	Mayo Stand Cover 80x145cm	1	EA



Laparoscopy Packs

NO.	BOM	Qty	Unit
1	Standard Surgical Gowns L	2	EA
2	Standard Surgical Gowns XL	2	EA
3	Laparoscopy Drape 320x210x350cm	1	EA
4	Hand Towels 30x40cm	4	EA
5	Mayo Stand Cover 80x145cm	1	EA
6	Back Table Cover 160x270cm	1	EA
7	OP-Tape 10x50cm	1	EA
8	Needle Counter 25 counts	1	EA
9	Tube Cover 15x250cm	1	EA



Cesarean Packs

NO.	BOM	Qty	Unit
1	Caesarean Drape 200x300cm	1	EA
2	Hand Towels 40x40cm	2	EA
3	Instrument Table cover 150x200cm	1	EA
4	Baby blanket 75x90CM	1	EA
5	Wrapper 100x100cm	1	EA
6	Umbilical Cord Clamp 5cm	1	EA
7	Reinforced Surgical Gowns L	2	EA
8	Mayo Stand cover 80x145cm	1	EA
9	Kidney Dish 700ml	1	EA
10	Maternity Pads	1	EA
11	Placenta Basin 1000ml	1	EA
12	Impervious Stockinette	1	EA
13	Elastic Bandage	1	EA
14	Syringe ball M	1	EA



Extremity Packs

NO.	BOM	Qty	Unit
1	Impervious Stockinette 22x75cm	1	EA
2	Liquid Collection Pouch 50x72cm	1	EA
3	Side Drape 75x90cm	1	EA
4	Split Drape 100x100cm	1	EA
5	Mayo Stand Cover 80x145cm	1	EA
6	Extremity Drape 200x300cm	1	EA
7	Instrument Table Cover 150x240cm	1	EA
8	Elastic Bandages 10x450cm	2	EA
9	Standard Surgical Gowns L	2	EA



Delivery Packs

NO.	BOM	Qty	Unit
1	Round Basin 8000ml	1	EA
2	Umbilical Cord Clamp 5cm	1	EA
3	Hand Towels 30x40cm	4	EA
4	Maternity Pad	1	EA
5	Leggings 75x120cm	2	EA
6	Elastic Bandages 10x450cm	2	EA
7	Reinforced Surgical Gowns L	2	EA
8	Instrument Table Cover 80x145cm	1	EA
9	Kidney Dish 700ml	1	EA
10	Mayo Stand Cover 80x145cm	1	EA
11	Suture Bag 25x30cm	1	EA
12	Bulb syringe M	1	EA
13	Abdominal Drape 230x270cm	1	EA
14	Baby Blankets 75x90cm	2	EA



Joint Procedure Packs

NO.	BOM	Qty	Unit
1	Instrument table covers 150x200 cm	2	EA
2	Mayo stand cover 80x145 cm	1	EA
3	Spilt Drape 200x260 cm	1	EA
4	Side drape 150x240 cm	1	EA
5	Side drape 150x200 cm	1	EA
6	Side drape 75x90 cm	1	EA
7	Leggings 40x120 cm	1	EA
8	Hand towels 40x40cm	4	EA
9	OP-tapes 10x50 cm	2	EA
10	Wrapper 100x100 cm	1	EA
11	Reinforced Surgical Gowns L	2	EA



Orthopedic Packs

NO.	BOM	Qty	Unit
1	Reinforced Surgical Gowns L	3	EA
2	Back Table Cover 150x200cm	1	EA
3	OP-Tape 10x50cm	1	EA
4	Side Drapes 150x220cm	2	EA
5	Hand Towels 30x40cm	3	EA
6	Instrument Table Cover 80x145cm	1	EA
7	Round Basin	1	EA
8	Needle Counter 25 counts	1	EA
9	U-Drape 200x290cm	1	EA
10	Liquid Collecting Pouch 35x45cm	1	EA



Angiography Packs

NO.	BOM	Qty	Unit
1	Kidney Dish 700ml	1	EA
2	Reinforced Surgical Gowns L	3	EA
3	Side Drapes 75x90cm	2	EA
4	Instrument Table Cover 80x145cm	1	EA
5	Basin XL	1	EA
6	Elastic Bandages 10x450cm	2	EA
7	Leggings 75x120cm	2	EA
8	Liquid Collecting Bag 35x45cm	1	EA
9	Mayo Stand Cover 150x200cm	1	EA
10	Medicine Cups 60ml	2	EA
11	OP Tape 10x50cm	1	EA
12	Cardiovascular Drape 232x380cm	1	EA



Arthroscopy Packs

NO.	BOM	Qty	Unit
1	Hand Towels 30x34cm	4	EA
2	Reinforced Surgical Gowns 125x150cm	2	EA
3	Arthroscopy Drape 200x320cm	1	EA
4	Liquid Collection Pouch 45x75cm	1	EA
5	Legging 40x75cm	1	EA
6	Skin Marker And Ruler 0.3-0.5cm	1	EA
7	OP Tapes 9x50cm	2	EA
8	Instrument Table Cover 150x200cm	1	EA
9	Wrapper 100x100cm	1	EA
10	Elastic Bandage 10x450cm	1	EA



TUR / Urology Packs

NO.	BOM	Qty	Unit
1	Instrument Table Cover 150x200 cm	1	EA
2	Urology Drape 180x240 cm	1	EA
3	Hand Towels 40x40 cm	2	EA
4	OP-Tapes 10x50 cm	2	EA
5	Wrapper 100x100 cm	1	EA
6	Standard Surgical Gowns L	2	EA



Dialysis Kits

NO.	BOM	Qty	Unit
1	Latex Examination Gloves XL	1	EA
2	Sterile Pad Sheet 50x50cm	1	EA
3	Tweezers Camel type/blue 13x4cm	2	EA
4	Tweezers Camel type/yellow 13x4cm	1	EA
5	Plaster 8.8x3cm	1	EA
6	Surgical Drape 60x60cm	1	EA
7	Shoe Cover 40x40cm	1	EA
8	Blue PVC tray	1	EA
9	Garbage Bag 25x30cm	1	EA
10	Cotton Balls 0.5g	5	EA
11	Tape 3.8x8cm	1	EA
12	Syringe 10ml	1	EA



Debridement kits

NO.	BOM	Qty	Unit
1	BandAids (1.5x2.3cm)	2	EA
2	Dry Cotton Balls	2	EA
3	Iodophor Cotton Ball	1	EA
4	PE Examination Gloves (1#)	1	EA
5	Tweezers	1	EA




Catheter Kits

NO.	BOM	Qty	Unit
1	Catheter	1	EA
2	Gauze Swab (5x7cm)	1	EA
3	Scrub	1	EA
4	Latex Examination Gloves	1	EA
5	Paraffin Cotton Ball	1	EA
6	Catheter Bowl	1	EA
7	Hole Drape	1	EA
8	Wrapper	1	EA
9	Drainage Bag	1	EA
10	Tray	1	EA
11	Tweezers	1	EA



Suture Kits

NO.	BOM	Qty	Unit
1	Latex Examination Gloves 7.5	1	EA
2	Apron 80x95cm	1	EA
3	Gauze Swab 5x5cm	5	EA
4	Hand Towel 30x40cm	1	EA
5	Sterile Pad Sheet 50x50cm	1	EA
6	Garbage Bag 25x30cm	1	EA
7	Cotton Balls 0.5g	5	EA
8	Tweezers Camel type/blue 13x4cm	2	EA
9	Blue PVC tray	1	EA
10	Syringe 10ml	1	EA
11	Paper Ruler 1.7x60cm	1	EA



Dressing Kits

NO.	BOM	Qty	Unit
1	Latex Examination Gloves XL	1	EA
2	Apron 80x95cm	1	EA
3	Gauze Swab 5x5cm	5	EA
4	Hand Towel 30x40cm	1	EA
5	Sterile Pad Sheet 50x50cm	1	EA
6	Garbage Bag 25x30cm	1	EA
7	Cotton Balls 0.5g	5	EA
8	Aseptic Dressing 10x10cm	1	EA
9	Tweezers Camel type/blue 13x4cm	2	EA
10	Blue PVC Tray	1	EA
11	Tweezers Camel type/yellow 13x4cm	1	EA



First aid Kits

NO.	BOM	Qty	Unit
1	Adhesive Bandages 3"x 1"	10	EA
2	Gauze Sponges 2"x 2"	4	EA
3	Alcohol Pads	4	EA
4	Povidone-iodine Pads	4	EA
5	Sting Relief Pads	2	EA
6	Cotton Tip Applicators	10	EA
7	Scissors	1	EA
8	Eye Wash 15ml	1	EA
9	Burn Gels 0.9g	2	EA
10	GKB104 PP Box	1	EA



I.V. Kits

NO.	BOM	Qty	Unit
1	Medical Tape 4.8cm	1	EA
2	Cotton Balls 0.5g	5	EA
3	Tourniquet 2x35cm	1	EA
4	Surgical Drape	1	EA
5	Hand Towel 30x40cm	1	EA
6	Blue PVC tray	1	EA
7	Tweezers Camel type/blue 13x4cm	2	EA
8	Tweezers Camel type/yellow 13x4cm	1	EA
9	Garbage Bag 25x30cm	1	EA
10	PVC inspection Glove XL	1	EA
11	Infusion Patch 42x80cm	1	EA



First aid Kits


NO.	BOM	Qty	Unit
1	PBT Bandage	1	EA
2	Knee/elbow Bandage	1	EA
3	Standard Adhesive Bandages	1	EA
4	Butterfly Bandages	1	EA
5	Medical Cotton Swab	1	EA
6	Emergency Foil Blankets	10	EA
7	Disposable CPR Masks	10	EA
8	Adhesive Wound Dressing	1	EA
9	Vinyl Gloves	1	EA
10	Sterile Eye Pads	4	EA
11	Manuals	4	EA
12	Alcohol Swabs	10	EA



Foley Catheters

Completely inert for less tissue irritation and encrustation during extended periods of indwelling use
 Easy insertion and less buckling with bullet-shaped tip
 Features symmetrical balloon shape for improved reliability and larger inner lumen for better drainage
 Clear for easy visualization of clots, mucous and urine flow
 Radiopaque strip provides easy visualization during scans

Silicone
 Fr8, Fr12, Fr14, Fr16, Fr18, Fr20, Fr22, Fr24
 Balloon 10ml, 30ml
 Not made with natural rubber latex



Connecting Catheters

Suction Connecting tubing is a clear, non-conductive tubing with two female ends used to connect aspirator to Yankauer suction device. Each tubing end has two soft universal female connectors for convenience and a secure attachment to any male connection.

1.75m, 2.45m, 2.95m, or Customized



Oxygen Tubes

The nasal cannula is a device used to deliver supplemental oxygen or increased airflow to a patient or person in need of respiratory help. The other end of the tube is connected to an oxygen supply such as a portable oxygen generator, or a wall connection in a hospital via a flowmeter.

Adult, Pediatric



Breathing Circuits

These breathing circuits are made of high quality tubing minimizing the risk of leakage for excellent durability. All connections are bonded together to prevent leakage and maintain the connection. All circuits are equipped with a secure hand grip allowing easy placement and removal from the anesthesia machine.

1.2m, 1.6m, 1.8m, or Customized



Vaginal Speculums

Smooth, molded PS design helps prevent tissue pinching or binding
Wider handle for improved maneuverability and balance
Strong grip and flexible yoke
Rack Fastener Type, Hook Fastener Type, Push-pull fastener Type

S M L

Not made with natural rubber latex



Gauze swabs

Produce and supply both with or without X-Ray detectable thread
Easy-open envelopes are carefully sealed to help block dust and contaminants
C-fold design helps minimize loose threads and lint
Cotton
7.5x7.5cm, 8x8cm, 10x10cm, or Customized



Surgical Gloves

Standard thickness, sterile, latex surgical gloves
Ideal for general surgical procedures
Made with natural rubber latex
Produce and supply both Powdered Gloves and Powder-Free Gloves
6.5, 7.0, 7.5, 8.0



Lap Sponges

Manufactured with absorbent cotton
Produce and supply both with or without handle and/or X-Ray detectable thread
Cotton
10x45cm, 30x30cm, 45x45cm, or Customized

Coveralls

Offer protection against infective agents (EN 14126), including resistance to penetration by blood and body fluids
Feature taped seams, an attached respirator fit hood (respirator fit), a storm flap over the zipper, and elastic wrists and ankles

Breathable Film

S M L XL XXL



Medical Face Masks

Type I, Type II, Type IIR, Level 1, Level 2
SPP+Meltblown+SPP
Blue, White or Customized
Earloop or Tie on
17.5x9.5cm



Medical Particulate Respirators

FFP1, FFP2, FFP3, KN95
SPP+Meltblown+Me.tbrown+SPP
Green, White or Customized



o laboratory/Workshop

